

Body Thermometer

You can modify the setting parameters to accommodate different features.

Accuracy: $\pm 0.2^{\circ}\text{C}$ / $\pm 0.4^{\circ}\text{F}$

- * Auto shut off 18 Seconds.
- * Measuring distance 1seconds.
- * Three-color backlight fever reminder.
- * You can store 50 sets of data.



Infrared Body Thermometer

Sample Rate: 0.5S
Accuracy: $\pm 0.2^{\circ}\text{C}$ (0.4°F)
Measuring: 1-5cm (0.4-2in)

Body
36.8

Healthy material shell

Use of healthy and environmentally friendly materials



ABS non-toxic environmentally friendly materials used to be assured



Ergonomic design
One-button press



Package



Packaging box




User's Manual



Body Thermometer

LCD Indicator

You can store 50 sets of data.



- | | |
|---------------------|---------------------------|
| ① Surface Mode | ② Body Mode |
| ③ Digital Display | ④ Smile or Cry Indicator |
| ⑤ Battery Indicator | ⑥ Buzzer ON/OFF Indicator |
| ⑦ Memory number | ⑧ °C ⑨ °F ⑩ Memories |



Infrared Probe

Measure the temperature by absorbing infrared heat emitted from the body

